
# SUSS MA6/8Gen4 光刻机

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General	Alignment	Alignment Views	Exposure
<b>Mask</b> Mask ID <input type="text" value="42"/> Mask Thickness [ $\mu\text{m}$ ] <input type="text" value="4700"/>		<b>Wafer</b> Wafer Thickness [ $\mu\text{m}$ ] <input type="text" value="725"/> Wafer Vacuum Minimum [mbar] <input type="text" value="-1000"/> Wafer Vacuum Maximum [mbar] <input type="text" value="-500"/>	
<b>Tool Set</b> Tool Set <input type="text" value="System 1 W150 M7 VAC"/> Chuck Vacuum Switch Sequence <input type="text" value="Inside"/>		<b>WEC</b> WEC Mode <input type="text" value="Wafer - Spacer - Top Substrate"/> WEC Offset [ $\mu\text{m}$ ] <input type="text" value="0"/>	


- Mask ID----Mask 名称
- Mask thickness----Mask 厚度
- Tool set----下拉选择 Tool set
- Chuck vacuum switch sequence----下拉选择 Chuck 真空方式
- Wafer thickness----Wafer 厚度
- Wafer thickness minimum----Wafer 厚度最小值
- Wafer thickness maximum----Wafer 厚度最大值
- WEC mode----下拉选择 WEC 方式
- WEC offset----WEC 补偿

General	Alignment	Alignment Views	Exposure
<b>Alignment Setup</b> Alignment Mode <input type="text" value="Manual"/> Do Alignment Check <input type="checkbox"/> Alignment Gap [ $\mu\text{m}$ ] <input type="text" value="50"/> Alignment Side <input type="text" value="Top Side"/> Single Alignment Microscope <input type="checkbox"/> Primary Illumination <input type="text" value="Top Coaxial"/> Secondary Illumination <input type="text" value="None"/> Top Lens <input type="text" value="10x"/>		<b>Align Positions</b> Left Mic Align Position X [ $\mu\text{m}$ ] <input type="text" value="-49933.125"/> Left Mic Align Position Y [ $\mu\text{m}$ ] <input type="text" value="52.5"/> Right Mic Align Position X [ $\mu\text{m}$ ] <input type="text" value="32000"/> Right Mic Align Position Y [ $\mu\text{m}$ ] <input type="text" value="0"/> Stage Align Position X [ $\mu\text{m}$ ] <input type="text" value="0"/> Stage Align Position Y [ $\mu\text{m}$ ] <input type="text" value="0"/> Stage Align Position $\Theta$ [ $^{\circ}$ ] <input type="text" value="0"/>	

- Alignment mode----对位方式
- Do alignment check----对位后检查
- Alignment gap----对位间隙
- Alignment side----对位方向
- Single alignment microscope----单显微镜对位
- Primary illumination----主照明
- Secondary illumination----次级照明
- Top lens----物镜类型
- Left mic align position X/Y----左显微镜位置 X/Y 轴
- Right mic align position X/Y----右显微镜位置 X/Y 轴
- Stage align position X/Y/ $\Theta$ ----载物台位置 X/Y/ $\Theta$  轴

General	Alignment	Alignment Views	Exposure
<b>Left Microscope</b>		<b>Right Microscope</b>	
<b>Mask</b>		<b>Mask</b>	
Primary Illu [%]	<input type="text" value="5"/>	Primary Illu [%]	<input type="text" value="5"/>
Secondary Illu [%]	<input type="text" value="0"/>	Secondary Illu [%]	<input type="text" value="0"/>
Top Focus [ $\mu\text{m}$ ]	<input type="text" value="25502.1875"/>	Top Focus [ $\mu\text{m}$ ]	<input type="text" value="25401.640625"/>
Bottom Focus [ $\mu\text{m}$ ]	<input type="text" value="0"/>	Bottom Focus [ $\mu\text{m}$ ]	<input type="text" value="0"/>
Gain [%]	<input type="text" value="0"/>	Gain [%]	<input type="text" value="0"/>
Contrast	<input type="text" value="0"/>	Contrast	<input type="text" value="0"/>
Integration Time [s]	<input type="text" value="0.035"/>	Integration Time [s]	<input type="text" value="0.035"/>
<b>Wafer in Alignment Gap</b>		<b>Wafer in Alignment Gap</b>	
Primary Illu [%]	<input type="text" value="5"/>	Primary Illu [%]	<input type="text" value="5"/>
Secondary Illu [%]	<input type="text" value="0"/>	Secondary Illu [%]	<input type="text" value="0"/>
Top Focus [ $\mu\text{m}$ ]	<input type="text" value="22944.140625"/>	Top Focus [ $\mu\text{m}$ ]	<input type="text" value="22676.210937"/>
Bottom Focus [ $\mu\text{m}$ ]	<input type="text" value="0"/>	Bottom Focus [ $\mu\text{m}$ ]	<input type="text" value="0"/>
Gain [%]	<input type="text" value="0"/>	Gain [%]	<input type="text" value="0"/>
Contrast	<input type="text" value="0"/>	Contrast	<input type="text" value="0"/>
Integration Time [s]	<input type="text" value="0.035"/>	Integration Time [s]	<input type="text" value="0.035"/>
<b>Wafer in Exposure Gap</b>		<b>Wafer in Exposure Gap</b>	
Primary Illu [%]	<input type="text" value="50"/>	Primary Illu [%]	<input type="text" value="50"/>
Secondary Illu [%]	<input type="text" value="0"/>	Secondary Illu [%]	<input type="text" value="0"/>
Top Focus [ $\mu\text{m}$ ]	<input type="text" value="20000"/>	Top Focus [ $\mu\text{m}$ ]	<input type="text" value="20000"/>
Bottom Focus [ $\mu\text{m}$ ]	<input type="text" value="0"/>	Bottom Focus [ $\mu\text{m}$ ]	<input type="text" value="0"/>
Gain [%]	<input type="text" value="50"/>	Gain [%]	<input type="text" value="50"/>
Contrast	<input type="text" value="50"/>	Contrast	<input type="text" value="50"/>
Integration Time [ms]	<input type="text" value="16"/>	Integration Time [ms]	<input type="text" value="16"/>

- Left Microscope----左显微镜
- Mask----掩膜版
- Wafer in alignment gap----晶圆在对位间隙
- Wafer in exposure gap----晶圆在曝光间隙
- Right Microscope----右显微镜
  - ◆ Primary illu----主照明
  - ◆ Secondary illu----次级照明
  - ◆ Top focus----上方焦距
  - ◆ Bottom focus----下方焦距
  - ◆ Gain----增益值
  - ◆ Contrast----对比度
  - ◆ Integration time----积分时间

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


The screenshot shows the SUSS Microtec M281AB Gen4 software interface. The top status bar includes fields for Recipe Name (SpaceredWEC HardConta), Process Type (Litho), Contact Mode (HardContact), TSA/BSA (Top Side), Alignment Gap (50 μm), Exposure Gap (0 μm), Current Gap (50 μm), User (Operator), and date/time (03.03.2017 09:21:29). The main window is divided into a central view area and a right-hand control panel. The central view shows a wafer with two crosshair alignment marks. The right panel has sections for General, Manual Alignment (Joystick Axes Positions: X: 0.8, Y: 1.2, Z: 2.1), Mask Images (Overlay, Grab), View and Focus (Mask View, Wafer View), and Save Recipe (Save, Save As). An Expose button is also present. At the bottom, there are tabs for Process, Recipes, and Service, and Back/Next navigation buttons.

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